

medifa
the medical factory

Made in Germany



medifa
the medical factory

Customized stainless steel solutions for

- Pathology
- Anatomy
- Forensic
- Veterinary
- Clinical equipment
- Stainless steel cabinets
- Project business
- Highest hygiene requirements
- High precision
- Strong technology expertise

**Project management, 3D SolidWorks engineering
and production**



DIN EN ISO 13485:2016
DIN EN ISO 9001:2015
DIN EN ISO 14001:2015

medifa metall und medizintechnik GmbH

Kehler Straße 31
76437 Rastatt
Germany
fon +49 (0)7222 7839-609
fax +49 (0)7222 7839-88
info.medizintechnik@medifa.com

we care.

**medifa customized medical
stainless steel solutions**



we care.

www.medifa.com



The medifa healthcare group unites competence concerning operating rooms under one roof.

medifa is an owner-run and internationally active complete solution provider in the medical technology.


The product and service portfolio covers all operating room requirements from modular wall, door, ceiling and cabinet systems (RooSy® Room Systems) for the turnkey operating room, as well as mobile operating tables, operating table accessories and operating lights, right up to examination chairs and medical furniture. In addition, one of the largest German OEM suppliers in the medical technology sector also belongs to medifa.

The medifa healthcare group employs approximately 350 employees at


three locations within Germany. The planning, engineering and production take place exclusively in Germany, and with over 90% added value in the group of companies, the group guarantees fast reaction times.

The slogan “we care” means for us that we consistently orientate ourselves according to the requirements of our customers. Customers all over the world can rely on top-class consultation, processing and quality from a single source, made in Germany according to certified national and international quality standards.

- Modular operating room systems: walls, doors, ceilings and cabinets (RooSy)
- Mobile operating tables and accessories
- Operating lights
- Examination and treatment chairs
- Medical furniture and trolleys
- Customized medical solutions



Solutions examples




Special construction operating table

- Height adjustable operating table
- Careful patient positioning on a removable, perforated gel surface
- Gentle procedure based on water vapour wetting of the body
- Effective removal of tissue remnants into the tub
- Control panel for disinfectant metering unit, water supply, electricity and medical gases
- Maximum body load: 300kg
- All inside edges with radius 15 mm for good cleaning



Autopsy and dissection tables

- Easy body handling
- With water, power supply and suction
- Mobile or stationary
- Dimensions as requested by the customer
- In combination with body trays and fresh air supply units with lightning
- Accessories like arm or head rest



Lab operating tables for animals

- Height adjustable operating table
- With water, power supply and extractions
- With tray holders and magnifier lights
- For small, mid-size and big animals

