



America

CERTIFICATE

No. U8 102069 0003 Rev. 00

Holder of Certificate: medifa GmbH & Co. KG
 Industriestr. 5
 57413 Finnentrop
 GERMANY

Production Facility(ies): 102069

Certification Mark:



Product:

Model(s):

Operating tables

medifa 8000

Parameters:

Voltage:	100 V - 240 V, 50 – 60 Hz
Battery:	24 VDC
Power:	300 VA
Protection class:	I
Safety working load:	250 kg
Weight:	250 kg
Applied part:	Type B
Duty cycle:	1/9 min.
Degree of protection:	IPX4

Tested according to:

ANSI/AAMI ES60601-1:2005/(R)2012
 CAN/CSA-C22.2 No. 60601-1:14

The product was voluntarily tested according to the relevant safety requirements noted above. It can be marked with the certification mark above. The mark must not be altered in any way. This product certification system operated by TÜV SÜD America Inc. most closely resembles system 3 as defined in ISO/IEC 17067. Certification is based on the TÜV SÜD "Testing and Certification Regulations". TÜV SÜD America Inc. is an OSHA recognized NRTL and a Standards Council of Canada accredited certification body.

Test report no.: 713123590

Date, 2019-01-10

(Stefan Reis)