

IEC SYSTEM FOR MUTUAL RECOGNITION OF TEST CERTIFICATES FOR ELECTRICAL EQUIPMENT
(IECEE) CB SCHEME

CB TEST CERTIFICATE

Product

Name and address of the applicant

Operating tables
medifa GmbH & Co. KG
Industriestr. 5
57413 Finnentrop
GERMANY

Name and address of the manufacturer

medifa GmbH & Co. KG
Industriestr. 5, 57413 Finnentrop, GERMANY

Name and address of the factory

medifa GmbH & Co. KG
Industriestr. 5, 57413 Finnentrop, GERMANY

Ratings and principal characteristics

Voltage:	100 V - 240 V, 50 – 60 Hz
Battery:	24 VDC
Power:	300 VA
Protection class:	I
Safety working load:	250 kg
Weight:	250 kg
Applied part:	Type B
Duty cycle:	1/9 min.
Degree of protection:	IPX4

Trade mark (if any)

medifa

Model/type Ref.

medifa 8000

A sample of the product was tested and found
to be in conformity with

IEC 60601-1:2005
IEC 60601-1:2005/AMD1:2012
IEC 60601-1-6:2010
IEC 60601-1-6:2010/AMD1:2013
IEC 60601-2-46:2016

as shown in the Test Report Ref. No.
which forms part of this certificate

713123590

This CB Test Certificate is issued by the National Certification Body

CB 102069 0002 Rev. 00

Date, 2019-01-08

(Stefan Reis)



Product Service

Page 1 of 1

TÜV SÜD Product Service GmbH • Certification Body • Ridlerstraße 65 • 80339 Munich • Germany